

Product Datasheet

LGX 818 (Encorafenib) (orb1183000)

Description	LGX 818 (Encorafenib)
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{22}H_{27}ClFN_7O_4S$
Purity	> 98%
MW	540.01
CAS Number	[1269440-17-6, 1269440-29-0]
SMILES	<chem>C[C@@H](CNC1=NC=CC(=N1)C2=CN(N=C2C3=CC(=CC(=C3F)NS(=O)(=O)C)Cl)C(C)C)NC(=O)OC</chem>
Expiration Date	12 months from date of receipt.