

Product Datasheet

GW6471 (orb1183299)

Description	GW6471
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{35}H_{36}F_3N_3O_4$
Purity	> 98%
MW	619.67
CAS Number	[436159-64-7]
SMILES	<chem>CCC(=O)NCC(CC1=CC=C(C=C1)OCCC2=C(OC(=N2)C3=CC=CC=C3)C)NC(=CC(=O)C4=CC=C(C=C4)C(F)(F)F)C</chem>
Expiration Date	12 months from date of receipt.