

Product Datasheet

LX 1606 (Telotristat) (orb1183425)

| | |
|------------------------|---|
| Description | LX 1606 (Telotristat) |
| Form/Appearance | Lyophilized |
| Storage | Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles. |
| Note | For research use only |
| Formula | $C_{27}H_{26}ClF_3N_6O_3$ |
| Purity | > 98% |
| MW | 575 |
| CAS Number | [1033805-22-9] |
| SMILES | <chem>CCOC(=O)[C@H](CC1=CC=C(C=C1)C2=CC(=NC(=N2)N)O[C@H](C3=C(C=C(C=C3)Cl)N4C=CC(=N4)C)C(F)(F)F)N</chem> |
| Expiration Date | 12 months from date of receipt. |