

Product Datasheet

Merck SIP Agonist (orb1183430)

Description	Merck SIP Agonist
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{23}H_{25}N_3O_3$
Purity	> 98%
MW	391.46
CAS Number	[635701-59-6]
SMILES	<chem>CC(C)CC1=CC=C(C=C1)C2=NC(=NO2)C3=CC=C(C=C3)CN4CC(C4)C(=O)O</chem>
Expiration Date	12 months from date of receipt.