

Product Datasheet

Tegaserod maleate (orb1183435)

Description	Tegaserod maleate
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{16}H_{23}N_5O \cdot C_4H_4O_4$
Purity	> 98%
MW	417.46
CAS Number	[189188-57-6]
SMILES	<chem>CCCCCN=C(N)NN/C=C/1C=NC2=C1C=C(C=C2)OC.C(=CC(=O)O)C(=O)O</chem>
Expiration Date	12 months from date of receipt.