

Product Datasheet

Lornoxicam (Xefo) (orb1183447)

Description	Lornoxicam (Xefo)
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{13}H_{10}ClN_3O_4S_2$
Purity	> 98%
MW	371.82
CAS Number	[70374-39-9]
SMILES	<chem>CN1C(=C(C2=C(S1(=O)=O)C=C(S2)Cl)O)C(=O)NC3=CC=CC=N3</chem>
Expiration Date	12 months from date of receipt.