

## Product Datasheet

### Phenylephrine HCl (orb1183457)

<b>Description</b>	Phenylephrine HCl
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>9</sub> H <sub>13</sub> NO <sub>2</sub> ·HCl
<b>Purity</b>	> 98%
<b>MW</b>	203.67
<b>CAS Number</b>	[61-76-7]
<b>SMILES</b>	<chem>CNC[C@@H](C1=CC(=CC=C1)O)O.Cl</chem>
<b>Expiration Date</b>	12 months from date of receipt.