

## Product Datasheet

### Rilpivirine (R 278474, TMC 278) (orb1183480)

<b>Description</b>	Rilpivirine (R 278474, TMC 278)
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{22}H_{18}N_6$
<b>Purity</b>	> 98%
<b>MW</b>	366.4
<b>CAS Number</b>	[500287-72-9]
<b>SMILES</b>	<chem>CC1=CC(=CC(=C1NC2=NC(=NC=C2)NC3=CC=C(C=C3)C#N)C)/C=C/C#N</chem>
<b>Expiration Date</b>	12 months from date of receipt.