

Product Datasheet

MLN4924 (HCL Salt) (orb1183488)

Description	MLN4924 (HCL Salt)
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{21}H_{26}ClN_5O_4S$
Purity	> 98%
MW	480
CAS Number	[1160295-21-5]
SMILES	<chem>C1CC2=CC=CC=C2[C@H]1NC3=NC=NC4=C3C=CN4[C@@H]5C[C@H]([C@H](C5)O)C(=O)(=O)N.Cl</chem>
Expiration Date	12 months from date of receipt.