

Product Datasheet

Glucocorticoid receptor agonist (orb1183504)

Description	Glucocorticoid receptor agonist
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{20}H_{20}F_4N_2O_2$
Purity	> 98%
MW	396.38
CAS Number	[1245526-82-2]
SMILES	<chem>OC1=CC=C(C=C1C(C)(C[C@](O)(C(F)(F)F)CC2=CC3=C(C=NC=C3)N2)C)F</chem>
Expiration Date	12 months from date of receipt.