

## Product Datasheet

### TAK-700 (Orteronel) (orb1183584)

<b>Description</b>	TAK-700 (Orteronel)
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{18}H_{17}N_3O_2$
<b>Purity</b>	> 98%
<b>MW</b>	307.4
<b>CAS Number</b>	[566939-85-3, 426219-18-3]
<b>SMILES</b>	<chem>CNC(=O)C1=CC2=C(C=C1)C=C(C=C2)[C@]3(CCN4C3=CN=C4)O</chem>
<b>Expiration Date</b>	12 months from date of receipt.