

Product Datasheet

TAK-875 (Fasiglifam) (orb1183590)

| | |
|------------------------|---|
| Description | TAK-875 (Fasiglifam) |
| Form/Appearance | Lyophilized |
| Storage | Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles. |
| Note | For research use only |
| Formula | C ₂₉ H ₃₂ O ₇ S |
| Purity | > 98% |
| MW | 524.6 |
| CAS Number | [1000413-72-8] |
| SMILES | <chem>CC1=CC(=CC(=C1C2=CC(=CC=C2)COC3=CC4=C(C=C3)[C@@H](CO4)CC(=O)O)C)OCCCS(=O)(=O)C</chem> |
| Expiration Date | 12 months from date of receipt. |