

Product Datasheet

XL184 free base (Cabozantinib) (orb1183593)

Description	XL184 free base (Cabozantinib)
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{28}H_{24}FN_3O_5$
Purity	> 98%
MW	501.5
CAS Number	[849217-68-1]
SMILES	<chem>COC1=CC2=C(C=CN=C2C=C1OC)OC3=CC=C(C=C3)NC(=O)C4(CC4)C(=O)NC5=CC=C(C=C5)F</chem>
Expiration Date	12 months from date of receipt.