

Product Datasheet

Roscovitine (Seliciclib) (orb1183611)

Description	Roscovitine (Seliciclib)
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{19}H_{26}N_6O$
Purity	> 98%
MW	354.5
CAS Number	[186692-46-6]
SMILES	<chem>CC[C@H](CO)NC1=NC2=C(C(=N1)NCC3=CC=CC=C3)N=CN2C(C)C</chem>
Expiration Date	12 months from date of receipt.