

Product Datasheet

Bibf1120 (Nintedanib) (orb1183654)

Description	Bibf1120 (Nintedanib)
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{31}H_{33}N_5O_4$
Purity	> 98%
MW	539.6
CAS Number	[928326-83-4, 656247-17-5]
SMILES	<chem>CN1CCN(CC1)CC(=O)N(C)C2=CC=C(C=C2)N/C(=C3/C4=C(C=C(C=C4)C(=O)OC)NC3=O)/C5=CC=CC=C5</chem>
Expiration Date	12 months from date of receipt.