

Product Datasheet

Aciclovir (Acyclovir) (orb1183666)

Description	Aciclovir (Acyclovir)
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_8H_{11}N_5O_3$
Purity	> 98%
MW	225.2
CAS Number	[59277-89-3]
SMILES	<chem>C1=NC2=C(N1COCCO)NC(=NC2=O)N</chem>
Expiration Date	12 months from date of receipt.