

Product Datasheet

CD1a antibody (APC) (orb1183976)

Description	Mouse monoclonal antibody CD1a conjugated to APC
Species/Host	Mouse
Conjugation	APC
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human CD1a
Target	CD1A
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG1, κ
Clonality	Monoclonal
MW	34 kDa
Uniprot ID	P06126
Dilution Range	Flow-Cyt=10ul/Test,IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.