

Product Datasheet

CD235a (GYPA) antibody (PE/Cy7) (orb1184016)

Description	Mouse monoclonal antibody CD235a (GYPA) conjugated to PE/Cy7
Species/Host	Mouse
Reactivity	Human
Conjugation	PE/Cy7
Tested Applications	FC
Target	GYPA
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG2a, k
Clonality	Monoclonal
MW	14 kDa
Uniprot ID	P02724
Dilution Range	Flow-Cyt=
Expiration Date	12 months from date of receipt.