

Product Datasheet

CD11b antibody (APC/Cy7) (orb1184035)

Description	Mouse monoclonal antibody CD11b conjugated to APC/Cy7
Species/Host	Mouse
Reactivity	Human
Conjugation	APC/Cy7
Tested Applications	FC
Target	ITGAM
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG1, κ
Clonality	Monoclonal
MW	125 kDa
Uniprot ID	P11215
Expiration Date	12 months from date of receipt.