

Product Datasheet

Human CD90/Thy-1 Mouse Monoclonal Antibody (FITC) (orb1184045)

Description	Human CD90/Thy-1 Mouse Monoclonal Antibody (FITC)
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Immunogen	HEL cells
Target	THY1
Preservatives	0.01M PBS (pH7.4), 0.5% rAlbumin, 0.02% Proclin300.
Form/Appearance	Liquid
Concentration	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Clone Number	2D9
Antibody Type	Primary Antibody
MW	25-35 kDa
Uniprot ID	P04216

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Dilution Range	Flow-Cyt=10ul/Test
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)