

Product Datasheet

CD273 (B7-DC, PD-L2) antibody (PE) (orb1184047)

Description	Mouse monoclonal antibody CD273 (B7-DC, PD-L2) conjugated to PE
Species/Host	Mouse
Conjugation	PE
Tested Applications	FC, ICC, IF
Target	PDCD1LG2
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	29 kDa
Uniprot ID	Q9BQ51
Dilution Range	Flow-Cyt=,IF=1:20-200,ICC=1:20-200
Expiration Date	12 months from date of receipt.