

## Product Datasheet

### CD13 antibody (FITC) (orb1184051)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Mouse monoclonal antibody CD13 conjugated to FITC   |
| <b>Species/Host</b>        | Mouse   |
| <b>Conjugation</b>         | FITC  |
| <b>Tested Applications</b> | FC  |
| <b>Target</b>              | ANPEP   |
| <b>Preservatives</b>       | 0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300           |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Storage</b>             | Store at 2-8°C. Protect from light. Avoid freezing. |
| <b>Note</b>                | For research use only                               |
| <b>Isotype</b>             | Mouse IgG2a, k                                      |
| <b>Clonality</b>           | Monoclonal  |
| <b>MW</b>                  | 109 kDa   |
| <b>Uniprot ID</b>          | <a href="#">P15144</a>                              |
| <b>Dilution Range</b>      | Flow-Cyt=   |
| <b>Expiration Date</b>     | 12 months from date of receipt.                     |