

Product Datasheet

CD19 antibody (APC/Cy7) (orb1184106)

Description	Rat monoclonal antibody CD19 conjugated to APC/Cy7
Species/Host	Rat
Conjugation	APC/Cy7
Tested Applications	FC
Target	CD19
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Rat IgG2a,k
Clonality	Monoclonal
MW	90 kDa
Uniprot ID	P15391
Expiration Date	12 months from date of receipt.