

Product Datasheet

CD117(c-Kit)/Pacific Blue antibody (orb1184120)

Description	Rat monoclonal antibody to CD117(c-Kit)/Pacific Blue
Species/Host	Rat
Conjugation	Pacific Blue
Tested Applications	FC
Target	KIT
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Rat IgG2b, k
Clonality	Monoclonal
MW	105 kDa
Uniprot ID	P10721
Expiration Date	12 months from date of receipt.