

Product Datasheet

NK1.1 antibody (FITC) (orb1184132)

Description	Mouse monoclonal antibody NK1.1 conjugated to FITC
Species/Host	Mouse
Conjugation	FITC
Tested Applications	FC
Target	Klrb1c
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG2a, κ
Clonality	Monoclonal
MW	25 kDa
Uniprot ID	P27814
Dilution Range	Flow-Cyt=
Expiration Date	12 months from date of receipt.