

Product Datasheet

CD123 antibody (PE/Cy7) (orb1184172)

| | |
|----------------------------|--|
| Description | Mouse monoclonal antibody CD123 conjugated to PE/Cy7 |
| Species/Host | Mouse |
| Conjugation | PE/Cy7 |
| Tested Applications | FC |
| Target | IL3RA |
| Preservatives | 0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300 |
| Form/Appearance | Liquid |
| Storage | Store at 2-8°C. Protect from light. Avoid freezing. |
| Note | For research use only |
| Isotype | Mouse IgG1, k |
| Clonality | Monoclonal |
| MW | 40 kDa |
| Uniprot ID | P26951 |
| Dilution Range | Flow-Cyt= |
| Expiration Date | 12 months from date of receipt. |