

Product Datasheet

CD123 antibody (PerCP/Cy5.5) (orb1184173)

Description	Mouse monoclonal antibody CD123 conjugated to PerCP/Cy5.5
Species/Host	Mouse
Conjugation	PerCP/Cy5.5
Tested Applications	FC
Target	IL3RA
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG1, k
Clonality	Monoclonal
MW	40 kDa
Uniprot ID	P26951
Dilution Range	Flow-Cyt=
Expiration Date	12 months from date of receipt.