

Product Datasheet

CD49b antibody (FITC) (orb1184224)

Description	Rat monoclonal antibody CD49b conjugated to FITC
Species/Host	Rat
Conjugation	FITC
Tested Applications	FC
Target	ITGA2
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Rat IgM, k
Clonality	Monoclonal
MW	127 kDa
Uniprot ID	P17301
Dilution Range	Flow-Cyt=
Expiration Date	12 months from date of receipt.