

Product Datasheet

CD49d antibody (FITC) (orb1184231)

Description	Rat monoclonal antibody CD49d conjugated to FITC
Species/Host	Rat
Conjugation	FITC
Tested Applications	FC
Target	ITGA4
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Rat IgG2b, k
Clonality	Monoclonal
MW	111 kDa
Uniprot ID	P13612
Dilution Range	Flow-Cyt=
Expiration Date	12 months from date of receipt.