

Product Datasheet

H-2Dd antibody (PE) (orb1184330)

Description	Mouse monoclonal antibody H-2Dd conjugated to PE
Species/Host	Mouse
Conjugation	PE
Tested Applications	FC
Target	H2-D1
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG2a, k
Clonality	Monoclonal
MW	41 kDa
Uniprot ID	P14426
Dilution Range	Flow-Cyt=
Expiration Date	12 months from date of receipt.