

Product Datasheet

H-2K(D)/H2-D1 antibody (PE) (orb1184332)

Description	Mouse monoclonal antibody H-2K(D)/H2-D1 conjugated to PE
Species/Host	Mouse
Conjugation	PE
Tested Applications	FC
Target	H2-K1
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG2a, k
Clonality	Monoclonal
MW	42 kDa
Uniprot ID	P04223
Dilution Range	Flow-Cyt=
Expiration Date	12 months from date of receipt.