

Product Datasheet

CD107b antibody (APC) (orb1184350)

Description	Rat monoclonal antibody CD107b conjugated to APC
Species/Host	Rat
Conjugation	APC
Tested Applications	FC
Target	LAMP2
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Rat IgG2a, k
Clonality	Monoclonal
MW	45 kDa
Uniprot ID	P13473
Dilution Range	Flow-Cyt=
Expiration Date	12 months from date of receipt.