

Product Datasheet

MPO antibody (FITC) (orb1184413)

Description	Mouse monoclonal antibody MPO conjugated to FITC
Species/Host	Mouse
Conjugation	FITC
Tested Applications	FC
Target	MPO
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Clonality	Monoclonal
MW	84 kDa
Uniprot ID	P05164
Dilution Range	Flow-Cyt=
Expiration Date	12 months from date of receipt.