

Product Datasheet

CD138 (Syndecan-1) antibody (APC) (orb1184423)

Description	Mouse monoclonal antibody CD138 (Syndecan-1) conjugated to APC
Species/Host	Mouse
Conjugation	APC
Tested Applications	FC
Target	SDC1
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG1, k
Clonality	Monoclonal
MW	34 kDa
Uniprot ID	P18827
Dilution Range	Flow-Cyt=
Expiration Date	12 months from date of receipt.