

## Product Datasheet

### CD196 (CCR6) antibody (PE/Cy7) (orb1184433)

<b>Description</b>	Mouse monoclonal antibody CD196 (CCR6) conjugated to PE/Cy7
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	PE/Cy7
<b>Tested Applications</b>	FC
<b>Target</b>	CCR6
<b>Preservatives</b>	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Store at 2-8°C. Protect from light. Avoid freezing.
<b>Note</b>	For research use only
<b>Isotype</b>	Mouse IgG2b, κ
<b>Clonality</b>	Monoclonal
<b>MW</b>	41 kDa
<b>Uniprot ID</b>	<a href="#">P51684</a>
<b>Dilution Range</b>	Flow-Cyt=
<b>Expiration Date</b>	12 months from date of receipt.