

Product Datasheet

Ly-6G/Ly-6C antibody (PE/Cy5) (orb1184469)

Description	Rat monoclonal antibody Ly-6G/Ly-6C conjugated to PE/Cy5
Species/Host	Rat
Conjugation	PE/Cy5
Tested Applications	FC
Target	Ly6g5c
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Rat IgG2b, k
Clonality	Monoclonal
MW	14 kDa
Uniprot ID	Q5SRR4
Expiration Date	12 months from date of receipt.