

## Product Datasheet

### Ly-6G/Ly-6C antibody (PE/Cy7) (orb1184470)

<b>Description</b>	Rat monoclonal antibody Ly-6G/Ly-6C conjugated to PE/Cy7
<b>Species/Host</b>	Rat
<b>Conjugation</b>	PE/Cy7
<b>Tested Applications</b>	FC
<b>Target</b>	Ly6g5c
<b>Preservatives</b>	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Store at 2-8°C. Protect from light. Avoid freezing.
<b>Note</b>	For research use only
<b>Isotype</b>	Rat IgG2b, k
<b>Clonality</b>	Monoclonal
<b>MW</b>	14 kDa
<b>Uniprot ID</b>	<a href="#">Q5SRR4</a>
<b>Expiration Date</b>	12 months from date of receipt.