

## Product Datasheet

### CD19 antibody (APC/Cy5.5) (orb1184513)

<b>Description</b>	Rat monoclonal antibody CD19 conjugated to APC/Cy5.5
<b>Species/Host</b>	Rat
<b>Conjugation</b>	APC/Cy5.5
<b>Tested Applications</b>	FC
<b>Target</b>	CD19
<b>Preservatives</b>	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Store at 2-8°C. Protect from light. Avoid freezing.
<b>Note</b>	For research use only
<b>Isotype</b>	Rat IgG2a,k
<b>Clonality</b>	Monoclonal
<b>MW</b>	59 kDa
<b>Uniprot ID</b>	<a href="#">P15391</a>
<b>Dilution Range</b>	Flow-Cyt=
<b>Expiration Date</b>	12 months from date of receipt.