

Product Datasheet

IL-17A antibody (PE) (orb1184519)

Description	Mouse monoclonal antibody IL-17A conjugated to PE
Species/Host	Mouse
Conjugation	PE
Tested Applications	FC
Target	IL17A
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG1, k
Clonality	Monoclonal
MW	15 kDa
Uniprot ID	Q62386
Dilution Range	Flow-Cyt=
Expiration Date	12 months from date of receipt.