

## Product Datasheet

### CD74 antibody (PE) (orb1184554)

<b>Description</b>	Mouse monoclonal antibody CD74 conjugated to PE
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	PE
<b>Tested Applications</b>	FC
<b>Target</b>	CD74
<b>Preservatives</b>	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Store at 2-8°C. Protect from light. Avoid freezing.
<b>Note</b>	For research use only
<b>Isotype</b>	Mouse IgG1, k
<b>Clonality</b>	Monoclonal
<b>MW</b>	34 kDa
<b>Uniprot ID</b>	<a href="#">P04233</a>
<b>Dilution Range</b>	Flow-Cyt=
<b>Expiration Date</b>	12 months from date of receipt.