

Product Datasheet

CD79a antibody (PE/Cy7) (orb1184557)

Description	Mouse monoclonal antibody CD79a conjugated to PE/Cy7
Species/Host	Mouse
Conjugation	PE/Cy7
Tested Applications	FC
Target	CD79A
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG1, k
Clonality	Monoclonal
MW	22 kDa
Uniprot ID	P11912
Dilution Range	Flow-Cyt=
Expiration Date	12 months from date of receipt.