

Product Datasheet

CD57 antibody (PE/Cy5) (orb1184600)

Description	Mouse monoclonal antibody CD57 conjugated to PE/Cy5
Species/Host	Mouse
Conjugation	PE/Cy5
Tested Applications	FC
Target	B3GAT1
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgM
Clonality	Monoclonal
MW	38 kDa
Uniprot ID	Q9P2W7
Dilution Range	Flow-Cyt=
Expiration Date	12 months from date of receipt.