



Product Datasheet

Annexin A10/ANXA10 Antibody (orb1184733)


Description	Annexin A10/ANXA10 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Immunogen	E.coli-derived human Annexin A10/ANXA10 recombinant protein (Position: N49-Y324).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	36 kDa
Uniprot ID	Q9UJ72
Dilution Range	Western blot, 0.25-0.5 µg/ml/ml, Human Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml/ml, Human Flow Cytometry, 1-3 µg/ml/1x10 ⁶ cells, Human Direct ELISA, 0.1-0.5 µg/ml/ml, Human
Expiration Date	12 months from date of receipt.



IHC analysis of Annexin A10/ANXA10 using...



IHC analysis of Annexin A10/ANXA10 using...



IHC analysis of Annexin A10/ANXA10 using...