



## Product Datasheet

### ME1 Antibody (monoclonal, 5E5F7) (orb1184749)


<b>Description</b>	ME1 Antibody (monoclonal, 5E5F7)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, WB
<b>Immunogen</b>	E.coli-derived human ME1 recombinant protein (Position: M1-Q572).
<b>Form/Appearance</b>	Lyophilized
<b>Concentration</b>	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
<b>Isotype</b>	Mouse IgG1
<b>Clonality</b>	Monoclonal
<b>MW</b>	64 kDa
<b>Uniprot ID</b>	<b>P48163</b>
<b>Expiration Date</b>	12 months from date of receipt.



IHC analysis of ME1 using anti-ME1 antib...



IHC analysis of ME1 using anti-ME1 antib...



IHC analysis of ME1 using anti-ME1 antib...