
Product Datasheet

Cynomolgus IL-2 R gamma / CD132 Protein (orb1184771)

Description

Cynomolgus IL-2 R gamma, His Tag is

(HEK293). It contains AA Leu 23 -...

Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-2 R gamma
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 29.4 kDa. The protein migrates as 50-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Cynomolgus IL-2 R gamma, His Tag (orb1184771) is expressed from human 293 cells (HEK293). It contains AA Leu 23 - Pro 255 (Accession # Q38JL2-1).
MW	29.4 kDa
Uniprot ID	Q38JL2
NCBI	NP_001030606.1
Expiration Date	6 months from date of receipt.