

Product Datasheet

Human BST1 Protein (orb1184784)

Description	Human BST1, His Tag is expressed from human 293 cells (HEK293). It contains AA Gly 29 - Ala 293 (Acc...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	BST1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 32.1 kDa. The protein migrates as 35-45 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_004325.2
Purity	95%
Source	Human BST1, His Tag (orb1184784) is expressed from human 293 cells (HEK293). It contains AA Gly 29 - Ala 293 (Accession # Q10588-1).
MW	32.1 kDa
Uniprot ID	Q10588
NCBI	NP_004325.2
Expiration Date	6 months from date of receipt.