

Product Datasheet

Human TRBC1 Protein (orb1184805)

Description	Human TRBC1, His Tag is expressed from human 293 cells (HEK293). It contains AA Asp 1 - Ala 144 (Acc...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	TRBC1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 18.2 kDa. The protein migrates as 20 kDa and 24 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Asp 1 - Ala 144
Purity	90%
Source	Human TRBC1, His Tag (orb1184805) is expressed from human 293 cells (HEK293). It contains AA Asp 1 - Ala 144 (Accession # P01850-1).
MW	18.2 kDa
Uniprot ID	P01850
Expiration Date	6 months from date of receipt.