

Product Datasheet

Rabbit Anti-Human IgG+IgM+IgA antibody (FITC) (orb1184934)

Description	Rabbit polyclonal antibody against Rabbit Anti-Human IgG+IgM+IgA conjugated to FITC
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Tested Applications	FC, IF
Immunogen	Human IgG+IgM+IgA
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Concentration	1mg/mL
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IF (1/250 - 1/2000), FC (1/250 - 1/2000)
Expiration Date	12 months from date of receipt.