

## Product Datasheet

# Rabbit Anti-Guinea Pig IgG (H&L) antibody (FITC) (orb1184944)

<b>Description</b>	Rabbit polyclonal antibody against Rabbit Anti-Guinea Pig IgG (H&L) conjugated to FITC
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Guinea pig
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	FC, IF
<b>Immunogen</b>	Guinea Pig IgG
<b>Preservatives</b>	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
<b>Concentration</b>	1mg/mL
<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	IF (1/250 - 1/2000), FC (1/250 - 1/2000)
<b>Expiration Date</b>	12 months from date of receipt.