

Product Datasheet

Rabbit Anti-Guinea Pig IgG (H&L) antibody (Biotin) (orb1184964)

Description	Rabbit polyclonal antibody against Rabbit Anti-Guinea Pig IgG (H&L) conjugated to Biotin
Species/Host	Rabbit
Reactivity	Guinea pig
Conjugation	Biotin
Tested Applications	ELISA, IHC, WB
Immunogen	Guinea Pig IgG
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Concentration	1mg/mL
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.